

**MEDICELL LTD**

A HAEMONETICS® COMPANY

# BASIC GUIDE TO TEG® INTERPRETATION



**MEDICELL LTD**  
DIAGNOSTICS DIVISION

**(t) 020 8371 9908**

# Cardiac Protocol

All samples are kaolin activated

Sample No.	When	Cup Type	Why Do It?
1	On induction	<ul style="list-style-type: none"> <li>Plain (clear) cup and pin</li> <li>If heparin suspected use heparinase bonded (blue) cup and pin as part of a split sample</li> </ul>	Gives baseline haemostatic profile
2	Re-Warmed (approx 36°C) (2)	<ul style="list-style-type: none"> <li>Heparinase bonded cup and pin (blue)</li> </ul>	Early identification of any coagulopathy that starts to develop during bypass
3	10 mins post protamine	Split sample: <ul style="list-style-type: none"> <li>Plain (clear) cup and pin</li> <li>Heparinase bonded cup and pin (blue)</li> </ul>	<ol style="list-style-type: none"> <li>Checks heparin reversal</li> <li>Gives post bypass haemostatic profile – identifies cause of coagulopathy in bleeding patient (see suggested treatment protocol)</li> </ol>
4	Post op ITU – if indicated	Split sample: <ul style="list-style-type: none"> <li>Plain (clear) cup and pin</li> <li>Heparinase bonded cup and pin (blue)</li> </ul>	<ol style="list-style-type: none"> <li>If it is indicated by a developing coagulopathy profile seen in samples 1–3</li> <li>If patient is bleeding or ‘going off’ clinically</li> </ol>

NB

- When using two cups ensure Plain cup is on the left of the heparinase cup
- This is the point the patient comes off the pump/ end of last anastomosis in off pump case.



(t) 020 8371 9908

# Suggested Treatment Protocol

## For kaolin activated samples

TEG® Values	Clinical Cause	Suggested Treatment
<b>R</b> between 11 – 14 min (2)	↓↓ clotting factors	X2 FFP or 8ml/kg
<b>R</b> greater than 14 min (2)	↓↓↓ clotting factors	X4 FFP or 16 ml/kg
<b>MA</b> between 42 – 47mm <b>G</b> between 3.6 – 4.4 K (3)	↓ platelet function	X 1 Platelet Pool
<b>MA</b> < 42 mm <b>G</b> less than 3.6K (3)	↓↓ platelet function	X 2 Platelet Pools
LY30 at 7.5% or greater, with CI less than 1.0	Primary fibrinolysis	Antifibrinolytic of choice
LY30 at 7.5% or greater, with CI greater than 3.0	Secondary fibrinolysis	Anticoagulant of choice
<b>R</b> less than 3 min <b>MA</b> greater than 75mm <b>G</b> greater than 15K	Prothrombotic state	Anticoagulant of choice


### NB

1. Blood product and antifibrinolytic treatment is only suggested if your patient is bleeding at a level you wish to treat
2. These guidelines apply after you have ruled out heparin as a cause for the long R time
3. A normal kaolin activated MA does not rule out platelet inhibitor effects. If the patient is/was on antiplatelet agents PlateletMapping™ results should be used to assess platelet function and hence need for treatment with platelets



(t) 020 8371 9908

# What Does TEG<sup>®</sup> Report?




# TEG<sup>®</sup> Parameter Definitions

- R = Reaction time
  - The reaction time is the time from when the blood was placed in the TEG<sup>®</sup> analyser until the initial fibrin formation ie how long it takes for the blood to start to clot, and is therefore a measure of clotting factors
- K = Kinetic time
  - The K time is a measure of the speed taken to reach a specific level of clot strength. Together with the alpha angle it is a measure of clot kinetics
- Alpha Angle
  - Measures the speed of fibrin build-up and cross-linking (clot strengthening)
- MA = Maximum Amplitude
  - The maximum amplitude represents the ultimate strength of the clot and is a measure of platelet function. It is a direct function of the maximum dynamic properties of fibrin and platelet bonding via GPIIb/IIIa.
- LY30 = Lysis 30 minutes after MA
  - It shows any breakdown of the clot and therefore gives an idea of clot stability




(t) 020 8371 9908

# TEG<sup>®</sup> Directed Therapy For Bleeding Patient




# Qualitative Analysis



# PlateletMapping™

There is a kit for monitoring the effect of platelet inhibiting drugs such as aspirin and clopidogrel.



The traces run for a maximum of 20 minutes. Overlaying the traces causes the software to calculate and display % inhibition and % function.



(t) 020 8371 9908



## What To Look For In An Etest

- Numbers in min/max are changing
- Final value 1800 – 2300
- Message etest ok
- If message is 'not at equilibrium' – retry
- If message is 'etest is off centre/out of range' – can be tweaked with trimmer tool – report to Site Administrator



(t) 020 8371 9908

# MEDICELL LTD

**70-72 Ballards Lane (Rear)  
The Mews, Off Princes Avenue, London N3 2BU  
T: +44 (0)20 8371 9908  
F: +44 (0)20 8371 9906**

**E: [enquiries@medi-cell.co.uk](mailto:enquiries@medi-cell.co.uk)  
[www.medi-cell.co.uk](http://www.medi-cell.co.uk)**

**MEDICELL LTD**

A HAEMONETICS® COMPANY

ME59011.REV04 – MediceLL Basic Guide to TEG® Interpretation



**(t) 020 8371 9908**